
Accurate Results for Textile Quality Control

Crock Meter

Manual Colour Fastness & Rubbing Fastness Tester

Digital Display | **9N** Weight | **Easy** Clamping | **Durable** Construction | **High** Accuracy



Digital Display | **9N** Weight | **Easy** Clamping | **Durable** Construction | **High** Accuracy

CROCK METER: Technical Specifications & Application

Features


- ✓ Digital Display
- ✓ 9N Standard Weight
- ✓ Easy Specimen Clamping
- ✓ Durable Structure
- ✓ Accurate Test Results
- ✓ Smooth Linear Slide Mechanism

Specifications

	N	N/mm	Cycle
Method	dry & wet rubbing method		
Weight	9N		
Operation	manual hand crank		
Counter	6-digit digital counter		
Compatible	ISO/AATCC compatible		

MODEL NO

GT-D05




Applications

- Textile Industry
- Garments Testing Laboratory Germents Department
- Quality Control Department
- Research Institutes
- University Laboratories

- Adopt Linear Slide, Running More Smooth.
- Rubbing Head Always Keep Horizontal, No Matter Thickness Of Specimen,No Adjust Position Of Rubbing Head, Easier Operation.
- LED Display Testing Times





LAB ASIA Science & Technology Corporation
THE RESULTS YOU CAN TRUST
www.labasiabd.com
labasia2020@gmail.com