LAB ASIA
Science & Technology Corporation
THE RESULTS YOU CAN TRUST

STAINLESS STEEL CALIBRATION WEIGHTS

High Quality | Accurate | Reliable



PRODUCT INFORMATION

Material	304 Stainless Steel
Class	F1, F2, M1
Nominal Value	1g,2g,5g,10g,20g,50g,100g,200g,500g,1kg, 2kg,5kg,10kg,20kg
Density	7.9g/cm^3
Magnetization	F1≤25μT; F2≤80μT; M1≤250μT
Package	Aluminum Box, Plastic Box
Susceptibility	0.05-0.2
Application extent	Apply to a high-precision laboratory, e.g.Metrology institution, Medicine factory, Food industry, and Environmental test institution.